BRAND	EQUIPMENT	APPROVED PRODUCT
1st Call Mobility	Mattresses	Clinell Universal
Abbott	FreeGo Pump	Clinell Universal
Abecca	Positioning Supports & Mattresses	Clinell Universal Clinell Detergent Clinell Peracetic Acid Wipes*
Active Key GmbH & Co. KG	All Active Key MedicalKey™ products, including: AK-C4112, AK-C4412 AK-C7012 & AK-C7412	Clinell Universal Clinell Detergent Clinell Peracetic Acid Wipes* Clinell Alcohol
AGFA	DX-D products	Clinell Universal Clinell Peracetic Acid Wipes*
Agile Medical	All products	Clinell Universal Clinell Detergent Clinell Peracetic Acid Wipes*
AKV International	Upholstery Vinyl	Clinell Universal
Albyn Medical	Hermes Urodynamic (do not use wipes on monitor)	Clinell Universal
Align Technology	iTero Element	Clinell Universal
All Modul	Trolley Range	Clinell Universal
Ambu	Ambu Display Units & Ambu aView 2 Advance	Clinell Universal Clinell Detergent
Arcomed	Pumps	Clinell Universal Clinell Detergent
ArjoHuntleigh	Mattresses Reliant IS2 Cover (this is a mattress cover)	Clinell Universal
ASep Healthcare	Tournistretch	Clinell Universal

\*Previously Clinell Sporicidal Wipes.

Approved Manufacturers. © GAMA Healthcare Ltd. Version 28.

www.gamahealthcare.com/compatibility

JBN230501



BRAND	EQUIPMENT	APPROVED PRODUCT
Baxter Healthcare	Evo IQ Infusion Pumps: ELVP001ANZ & ESYRANZ001	Clinell Universal
	Accessories: EHAN010, EHUBB005ANZ, EHUB006ANZ, EHUB007ANZ & EHUB008ANZ	
BD (Bard)	CareFusion Jaeger Spirometry Cabinet Alaris Syringe Pump Alaris GW Volumetric Pump	Clinell Universal
	Site~Rite® 8 Ultrasound System and Probe	Clinell Universal Clinell Peracetic Acid Wipes*
Bedfont	Smokerlyzer, ToxCO (not suitable with Gastrolyzer)	Clinell Universal
Brandon Medical	Astralite Treatment Lights Astramax Minor Surgical Lights Quasar Elite Operating Theatre Lights	Clinell Universal
Braun	Thermoscans	Clinell Universal
Carrflex	Mattress Covers, Upholstery Tops & Associated Equipment	Clinell Universal Clinell Detergent Clinell Peracetic Acid Wipes*
Carr Plastics	Techmaflex	Clinell Universal Clinell Detergent Clinell Peracetic Acid Wipes*
Central Medical	Swabsafe Stands	Clinell Universal
Covidien	Genius 2 Tympanic Thermometer	Clinell Universal Clinell Detergent

\*Previously Clinell Sporicidal Wipes.

Approved Manufacturers. © GAMA Healthcare Ltd. Version 28.

www.gamahealthcare.com/compatibility

JBN230501



BRAND	EQUIPMENT	APPROVED PRODUCT
Dartex Coating	Mattresses (with polyurethane face)	Clinell Universal Clinell Peracetic Acid Wipes*
De Smit Medical	CubeScan Bladder Scanners	Clinell Universal
Direct Healthcare	Dartex Mattress Covers	Clinell Universal
Echosens	FibroScan	Clinell Universal
Elekta	IBeam Evo Couch Top	Clinell Universal
Fresenius Kabi - Agilia range	VP MC Volumetric Infusion Pump	Clinell Universal
	SP MC Syringe Infusion Pump	Clinell Detergent
	SP TIVA TCI Syringe Infusion Pump	Clinell Peracetic Acid Wipes*
	SP PCA Syringe Infusion Pump	
Fritz Stephen	Pediatric Ventilation Systems	Clinell Universal Clinell Detergent
Fujifilm Healthcare	Ultrasound Probes and Accessories (suitable for most, please see cleaning instruction QAD-Z3510 for full details)	Clinell Universal Clinell Peracetic Acid Wipes*
Fujifilm Sonosite	Ultrasound Equipment (Edge, NanoMaxx, MicroMaxx, TITAN) iViz Devices and C60v, L25v, L38v & P21v probe	Clinell Universal Clinell Peracetic Acid Wipes*

\*Previously Clinell Sporicidal Wipes.

Approved Manufacturers. © GAMA Healthcare Ltd. Version 28.

www.gamahealthcare.com/compatibility

JBN230501



BRAND	EQUIPMENT	APPROVED PRODUCT
GE Healthcare	Ultrasound Transducers (please email info@gamahealthcare.com for full list of compatible transducers) Probe Holders, User interface, Touch panel, OLED Monitor Display and the Housing of Voluson Expert BT15, BT16, BT17, BT18 and BT19 Consoles. Note: LCD Monitor Display can only be cleaned with IPA (70%) solution Voluson Ultrasound Consoles Expert Series (E6, E8 and E10) Signature Series (S8 and S10) Performance Series (P8) MAC VU360 <sup>™</sup> Resting ECG	Clinell Universal
Haines Medica Australia	Haines MediCurtains (full range) Haines Wipeclean Waterproof Pillows Haines Reusable Tourniquets Haines PVC Wipedown Mattress Covers Haines PVC Mackintosh Rolls Haines Reusable Pillow Cases	Clinell Universal
Hospedia	Patient Entertainment Unit	Clinell Universal
Hospira	Sapphire Pump	Clinell Universal
Howard Wright	Mattresses	Clinell Universal Clinell Detergent Clinell Peracetic Acid Wipes*

\*Previously Clinell Sporicidal Wipes.

Approved Manufacturers. © GAMA Healthcare Ltd. Version 28.

www.gamahealthcare.com/compatibility



BRAND	EQUIPMENT	APPROVED PRODUCT
HP	IPAQ HP EliteBook 840 G6	Clinell Universal Clinell Universal
	Healthcare Edition HP EliteOne 800 G5 Healthcare Edition All-in-One	Clinell Peracetic Acid Wipes*
	HP Healthcare Edition HC270cr Clinical Review Display	
	HP Healthcare Edition HC241 Clinical Review Monitor	
	HP Healthcare Edition HC241p Clinical Review Monitor	
	HP Healthcare Edition HC271 Clinical Review Monitor	
	HP Healthcare Edition HC271p Clinical Review Monitor	
Huntleigh Healthcare	Foetal Monitors	Clinell Universal Clinell Detergent
International Light Technologies	IL1350	Clinell Universal
Karomed	Mattresses	Clinell Universal
Konica Minolta	SONOVISTA FX Ultrasound Device (premium edition)	Clinell Universal
Laerdal	Suction Units	Clinell Universal Clinell Detergent
MePACS	Duress Alarms	Clinell Universal
Mindray	N Series Monitor	Clinell Universal Clinell Peracetic Acid Wipes*

\*Previously Clinell Sporicidal Wipes.

Approved Manufacturers. © GAMA Healthcare Ltd. Version 28.

www.gamahealthcare.com/compatibility

JBN230501



BRAND	EQUIPMENT	APPROVED PRODUCT
Modsel	All Product Ranges	Clinell Universal Clinell Detergent Clinell Peracetic Acid Wipes*
OES Medical	Pressure Monitor and UAM Anaesthetic Machine	Clinell Universal Clinell Detergent Clinell Peracetic Acid Wipes*
Orfit Industries	All Related Types of AIO 3.0	Clinell Universal Clinell Detergent Clinell Peracetic Acid Wipes*
	Medium Density Blue Head Support	Clinell Universal Clinell Detergent Clinell Peracetic Acid Wipes*
	Low Density Blue Head Support	Clinell Detergent Clinell Peracetic Acid Wipes*
	All Related Types of SBRT Cushions Clinell Detergent	
OxyLitre	All Equipment	Clinell Universal Clinell Detergent
Panaz	Cadet Furniture	Clinell Universal Clinell Peracetic Acid Wipes*
Parity Medical	Mobile Computer Carts	Clinell Universal Clinell Peracetic Acid Wipes*
Park House	Healthcare Mattress Covers	Clinell Universal
Phantom Laboratory	Catphan 500	Clinell Universal

\*Previously Clinell Sporicidal Wipes.

Approved Manufacturers. © GAMA Healthcare Ltd. Version 28.

www.gamahealthcare.com/compatibility

JBN230501



BRAND	EQUIPMENT	APPROVED PRODUCT
Philips	Ultrasound Equipment Intellivue Monitor (remaining residue must be wiped off) Ultrasound Transducers (please email info@gamahealthcare.com	Clinell Detergent Clinell Universal Clinell Universal Clinell Peracetic Acid Wipes*
	for full list of compatible transducers)	·
Plinth Medical	Vinyl Couches	Clinell Universal
Recticel	Mattresses	Clinell Universal Clinell Peracetic Acid Wipes*
Seca	All Equipment	Clinell Universal
Sepura	Sepura TETRA Products	Clinell Universal
Sidhil	IQ Bed Frame	Clinell Universal Clinell Peracetic Acid Wipes*
Siemens	Mammomat Inspiration Siemens RapidPoint Systems	Clinell Universal Clinell Peracetic Acid Wipes*
Smiths Medical	Pressure Infuser	Clinell Universal
Spes Medica	LTM autoclavable cup electrode Ag/AgCl LTM autoclavable cup electrode Gold Ear clip electrode Ag/AgCl Bridge electrodes Sintered Ag/AgCl Moulded cup electrodes Ag/AgCl Paediatric cup electrodes Gold Cup electrodes Sintered Ag/AgCl Reusable airflow thermal sensor Tremor sensor Limb movement sensor Reusable Snoring Sensor	Clinell Universal

\*Previously Clinell Sporicidal Wipes.

Approved Manufacturers. © GAMA Healthcare Ltd. Version 28.

www.gamahealthcare.com/compatibility

JBN230501



BRAND	EQUIPMENT	APPROVED PRODUCT
Toshiba	Ultrasound Equipment	Clinell Universal Clinell Peracetic Acid Wipes*
Tough-PAC	iPad Covers	All Clinell products
Triton Electronics Systems Ltd	Multigas Analyser-AMG-06	Clinell Universal Clinell Detergent
Universal Ferrule Remover	Ferrule Remover	Clinell Universal Clinell Peracetic Acid Wipes*
Verathon	GlideScope Cables BladderScan Prime	Clinell Universal
Visionflex	Thinklabs One Digital Stethoscope	Clinell Universal
Welch Allyn	Connex Spot Monitor Caretemp Thermometers Pro6000 Thermometers Pocket Plus LED Diagnostic Set - Ophthalmoscope	Clinell Universal
Zebra	Zebra Healthcare Mobile Computers	Clinell Universal
Zoll	AED	Clinell Universal Clinell Peracetic Acid Wipes*

\*Previously Clinell Sporicidal Wipes.

Approved Manufacturers. © GAMA Healthcare Ltd. Version 28.

www.gamahealthcare.com/compatibility

JBN230501

